

Comfort - Functionality - Durability

Features:

- The Lumex FR587W Recliner features three positions for patient comfort. The fourth position is a caregiver-activated Trendelenburg (or Shock) Position
- Improved wider and softer armrest for patient comfort
- Wide range of fabrics / colors
- Quick Ship available
- Dual gas cylinders give smooth power assistance during Trendelenburg activation
- Wide, ergonomically designed padded armrests for improved patient comfort
- Four Tente[®] swivel casters facilitate transport and maneuverability in confined areas. Three total locking casters with one directional locking caster

- Durable, high-impact side tables sit flush with recliner arm to conveniently hold prep materials. Side table surface is smooth and features tapered edges for easy cleaning
- 450 lb maximum weight capacity, evenly distributed

Options:

- Heat and Massage
- Fully Upholstered Armrests
- IV Pole Holder / IV Pole
- Headrest Cover
- Universal Pillow
- Drop Arm Model Available (FR587WD)

Dimensions

Overall Height.....	50" / 127 cm
Overall Width.....	33" / 84 cm
Overall Depth.....	36" / 91 cm
Width Between Arms.....	23" / cm
Arm Height From Seat.....	8" / 20 cm
Seat Height.....	21.5" / 55 cm
Seat Width.....	23.5" / 58 cm
Seat Depth.....	20.5" / 52 cm
Depth Fully Reclined In Trendelenburg.....	80" / 203 cm
Width With Both Side Tables in Upright Position.....	53.25" / 135cm
Maximum Weight Capacity Evenly Distributed.....	450 lb / 204.1 kg



Multiple Positions for Enhanced Patient Comfort